

Fully autoclaveable metal blade, covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.

The non-replaceable glass fibre bundle is fully enclosed in the blade which allows easy cleaning.

Rigid metal construction from fully autoclaveable, antimagnetic stainless steel.

Lever positioned at the rear of the blade to control the tilting tip.

Clear product information etched on the blade.

Atraumatic distal end.

Tilting tip for elevation of the epiglottis, providing an improved view of the larynx, without excessive force, especially for difficult intubations.


Green ISO 7376 fitting metal block

3 metal ball bearings for secure fitting.



#### DIFFICULT INTUBATIONS


The Sirius XL Eclipse laryngoscope blades are an invaluable aid for difficult intubations as they feature a flexible tip for elevation of the epiglottis fold providing an improved view of the trachea without excessive force resulting in less patient trauma.



#### SIRIUS XL ECLIPSE BLADES

Description	Measurements
Sirius XL Eclipse Macintosh Size 2	115mm x 24mm
Sirius XL Eclipse Macintosh Size 3	135mm x 24mm
Sirius XL Eclipse Macintosh Size 4	160mm x 24mm

#### INTUBATION USING A TIMESCO SIRIUS XL ECLIPSE BLADE



The Timesco Sirius XL Eclipse laryngoscope blades benefit from an articulating tip to enable the elevation of the epiglottis fold; providing an improved view of the trachea without excessive force, resulting in less risk of patient trauma especially in difficult intubations. For example:

- Patients with restricted mouth opening.
- Patients with protruding teeth.
- Patients that have decreased neck movement.
- Patients with Spinal injuries.
- Patients with cervical neck collars in use.

Timesco also offer single use Eclipse blades in our Callisto range. Callisto Eclipse brochure: [L.L.069/1](#)